



Digital Veterinary Thermometer for Big Animal



Weiyou® is a professional digital veterinary thermometer manufacturer and supplier in China. We have been specialized in veterinary products for many years. digital veterinary thermometer for big animals have competitive price advantage and cover most of the European and American markets. Weiyou® hope to cooperate with more customers for mutual development and benefits.

Product Description

Weiyou® digital veterinary thermometer for big animals featured at its high accuracy and safe measurement. Quick and accurate reading in 10 seconds facilitates the measurement for the veterinarian, vet technician or the owner. Series of digital veterinary thermometer for big animals are in compliance with several certificates such as CE, ISO13485, GMP, FCC, etc. Weiyou® digital veterinary thermometer for big animals is mainly used for Cattle, Horse Sheep, Goat and other animals.

Specification

Item	Digital veterinary thermometer for big animals
Size	230 x 30 x 16mm
Material	ABS
MOQ	100pcs
Packing	1 pc in one box

Features

Veterinary thermometer

This product is degree C and F

Battery included: one LR44 alkaline button cell

Probe: stainless steel, length 12.7cm, diameter 0.5cm, blunt tip

















Product Description

1. Measuring range: 32°C to 42°C
2. Accuracy: +/- 0.2 C°
3. Probe: Stainless steel, length 12.7cm, diameter 0.5cm, blunt tip
4. Battery: One LR44 alkaline button cell
5. Auto-shutoff after 8 minutes.
6. Provide lower voltage detector function
7. Provide fever alarm function



Digital Veterinary Thermometer

Fast & Accurate
Versatile & Safe





